

pNeuton mini

Pneumatic Neonatal / Infant Ventilation

Operational Verification Procedure pNeuton mini Ventilator

Thank you for your purchase. This procedure outlines a validation testing method used to verify proper operation of the **pNeuton mini** Ventilator. It must be performed prior to using the ventilator for the first time in a health care facility. If the ventilator fails to pass any of the following tests do not apply it to patients. Call Airon Corporation Customer Support at 888-448-1238.

Ventilator Setup

The following equipment is needed:

1. **pNeuton mini** Ventilator with breathing circuit (Airon # 58031, 58035 or 58201)
2. Oxygen analyzer
3. Watch, a timer or stop watch

When ready:

1. Attach breathing circuit to ventilator following instructions in the Operators Manual.
2. Set the controls as follows:
 - a. **Mode** control to IMV + CPAP
 - b. **PEEP / CPAP** to off
 - c. **Peak Pressure (P_{INSP})** to 30 cm H₂O
 - d. **Continuous Flow (V_{INSP})** to 8 L/min
 - e. **Inspiratory Time (T_I)** to 0.6 sec
 - f. **Expiratory Time (T_E)** to 2.0 sec
 - g. **Oxygen (F_IO₂)** to 60%
 - h. **High Pressure Alarm (P_{ALARM})** to 40 cm H₂O
3. Occlude the patient connection or place a cap over the patient connection to seal the circuit.
4. Attach Oxygen and Air Inputs on side panel of the ventilator to high-pressure gas sources (55 psi ± 15 psi each gas) and turn on both gases.

NOTE: The ventilator will begin operation at the above settings when the oxygen and air are turned on. The alarm will sound. You may press the "Reset / Silence" button to silence the alarm.

Operational Verification

Verification Step	Acceptable Range	Result
Observe the pressure gauge. It should rise to 30 cm H ₂ O and drop to 0 cm H ₂ O.	30 ± 3 cm H ₂ O	Pass / Fail
Count the respiratory rate with a stopwatch. Measure the number of breaths in one minute.	22 ± 5 breaths per minute	Pass / Fail
Attach an oxygen analyzer to the output of the expiratory valve. Using a calibrated oxygen analyzer measure the oxygen percentage.	60% ± 3%	Pass / Fail
Set the High Pressure Alarm to 20 cm H ₂ O. The alarm should activate with each breath. Turn High Pressure Alarm back to 40 cm H ₂ O post test.	Visual and audible alarm	Pass / Fail
Remove the patient connection cap and allow the breathing circuit to remain open. Using a stopwatch, measure the time until the alarm sounds.	10 ± 2 seconds	Pass / Fail
Disconnect one gas input with the remaining gas source connected to Air or Oxygen gas source. Verify the Low Gas Supply alarm activates.	Visual and audible alarm	Pass / Fail

If the ventilator has passed all the above steps, it is ready for clinical use.

Tested by: _____ Date: _____



Airon

Airon Corporation
751 North Drive, Unit 6, Melbourne, Florida 32934 USA
tel (321) 821-9433 fax (321) 821-9443
CD-I-003 Rev E