

# pNeuton mini NEO

## Pneumatic Neonatal Ventilation

### Operational Verification Procedure pNeuton mini NEO Ventilator

This procedure outlines a validation testing method used to verify proper operation of the **pNeuton mini NEO** Ventilator. The Operational Verification should be performed prior to using the ventilator with each patient application.

#### Ventilator Setup

The following equipment is needed:

1. **pNeuton mini neo** Ventilator with breathing circuit (Airon # 58031 or 58201)
2. Oxygen analyzer
3. Watch, a timer or stop watch

When ready:

1. Attach breathing circuit to ventilator following instructions in the Operators Manual.
2. Set the controls as follows:
  - a. **Mode** control to IMV + CPAP
  - b. **PEEP / CPAP** to off
  - c. **Peak Pressure (P<sub>INSP</sub>)** to 20 cm H<sub>2</sub>O
  - d. **Continuous Flow (V<sub>INSP</sub>)** to 8 L/min
  - e. **Inspiratory Time (T<sub>I</sub>)** to 0.4 sec
  - f. **Expiratory Time (T<sub>E</sub>)** to 2.0 sec
  - g. **Oxygen (F<sub>I</sub>O<sub>2</sub>)** to 60%
  - h. **High Pressure Alarm (P<sub>ALARM</sub>)** to 30 cm H<sub>2</sub>O
3. Occlude the patient connection or place a cap over the patient connection to seal the circuit.
4. Attach Oxygen and Air Input on side panel of the ventilator to high-pressure gas sources (55 psi ± 15 psi each gas) and turn on both gases.

NOTE: The ventilator will begin operation at the above settings when the oxygen and air are turned on. The alarm will sound. You may press the "Reset / Silence" button to silence the alarm.

#### Operational Verification

Verification Step	Acceptable Range	Result
Observe the pressure gauge. It should rise to 20 cm H <sub>2</sub> O and drop to 0 cm H <sub>2</sub> O.	20 ± 3 cm H <sub>2</sub> O	Pass / Fail
Count the respiratory rate with a stopwatch. Measure the number of breaths in one minute.	25 ± 5 breaths per minute	Pass / Fail
Attach an oxygen analyzer to the output of the expiratory valve. Using a calibrated oxygen analyzer measure the oxygen percentage.	60% ± 3%	Pass / Fail
Set the High Pressure Alarm to 20 cm H <sub>2</sub> O. The alarm should activate with each breath. Turn High Pressure Alarm back to 30 cm H <sub>2</sub> O post test.	Visual and audible alarm	Pass / Fail
Remove the patient connection cap and allow the breathing circuit to remain open. Using a stopwatch, measure the time until the alarm sounds.	10 ± 2 seconds	Pass / Fail
Disconnect one gas input with the remaining gas source connected to Air or Oxygen gas source. Verify the Low Gas Supply alarm activates.	Visual and audible alarm	Pass / Fail

If the ventilator has passed all the above steps, it is ready for clinical use.

Tested by: \_\_\_\_\_ Date: \_\_\_\_\_



**Airon**

Airon Corporation  
751 North Drive, Unit 6, Melbourne, Florida 32934 USA  
tel 321-821-9433 fax 321-821-9443  
CD-IN-008 Rev A