


Easy-to-Use Pediatric Resuscitation Guide

Most of us don't deal with pediatric patients that often. When we do, especially if it's a critically ill or injured kid, the stress level goes up, making it hard to recall age- and weight-specific treatments and dosages. The **Pediatric Resuscitation Guide** from **Huston-Haines** was designed by paramedics to help take the stress out of the decision-making process and give you clear and concise length-based guidelines to follow when treating a pediatric patient. Everything you need is on one page—you don't have to flip from front to back to get the information you're looking for. The guide also features expected 50th percentile vital signs for each size child and also lists size recommendations for defibrillation pads, intraosseous needles and endotracheal tubes. Guides can be customized with your service's logo and set to your service's protocols.

VITALS

Sizes: 8.5" x 11" (standard size);
7" x 5" (travel size)
Price: \$29.95
www.pediatricresuscitation.com
920-680-0904



Dual Blade Laryngoscope

When you need to perform an endotracheal intubation in the field, it's usually a matter of life and death. You only have seconds to decide what size and shape blade to use, find that size in your airway kit, assemble the equipment and place the tube. The **DUOSCOPE** from **Airon** is a compact, easy-to-use, two-blades-in-one disposable emergency laryngoscope. The DUOSCOPE requires no assembly—simply pick the size you want to use, open and turn on the light. The DUOSCOPE is available in three different models: a child model with a Macintosh 3 and Miller 1 (RB2000), an adult model with a Macintosh 4 and Miller 1 (RB3000), and a large adult model with Macintosh 4 and Macintosh 3 (RB4000). The DUOSCOPE is constructed of high-quality lightweight plastic and the dual illumination lamps allow for over 24 hours of continuous illumination.

VITALS

Aperture opening: 90 degrees
Batteries: Two 3-volt lithium coin cell (CR2032)
Price: Call for price
www.aironusa.com
888-448-1238



Keep Your Legs Dry

Although most of the EMS jackets available on the market do a great job of keeping the upper part of your body dry in nasty weather, it can be difficult to keep your legs dry. The **911 Rain Pants** from **Gerber Outerwear** are designed to not only keep you dry in bad weather, but they're also NFPA-certified bloodborne pathogen impervious. 911 Rain Pants are lightweight, breathable and windproof. They feature large cargo pockets on each thigh, with striping available in silver Scotchlite or lime/silver triple trim. The zippers on each cuff make the pants fast and easy to put on or take off, and the elastic drawstring waist with adjustable barrel clips ensures the pants stay in place.

VITALS

Color: Midnight navy
Price: \$240.00
www.gerberouterwear.com
800-GERBER1

IN THE NEXT ISSUE: >> EPIC-id USB Emergency ID Braclet >> MegaMover Select >>> ShockLink
>>> Mustang Traction Splint >>> Active Shooter Event Causality Response Kit >>> Simulated IV Bags