



FOR IMMEDIATE RELEASE
December 10, 2015

For more Information:
John Fitch, Marketing
jfitch@AironUSA.com
321-821-9433

Duoscope and Onescope Released for USA Distribution

Airon Corporation announces the release of the **Duoscope and Onescope Disposable Laryngoscopes**, the next generation laryngoscopes for EMS, military operations and multiple hospital settings including: emergency department, post anesthesia, and intensive care. The **Duoscope** is an innovative double blade laryngoscope with 2 different size intubation blades. The **Onescope** is a standard single blade laryngoscope.



The **Duoscope** is truly unique. It's functional design facilitates immediate accessibility and instantaneous response in emergency situations. When opened the user has 2 different intubation blades to use. One blade functions as the handle while the other blade is inserted for cord visualization. The **Duoscope** comes in three different configurations, providing a combination of 4 different blade sizes.

The **Onescope** looks like a standard laryngoscope with a single blade. However, the handle is hollow making it a very lightweight and easy to use laryngoscope.

Both laryngoscopes use a Super Bright Pure White Water Clear LED that will operate for up to 24 hours. The laryngoscopes meet the ISO standard for blade structural integrity: tested with an applied load of 40 kg force, there was no breakage for the full range of the blade.

The **Duoscope and Onescope Disposable Laryngoscopes** are available for purchase through Airon's domestic distribution partners or directly from Airon Corporation. For more information call 888-448-1238 (USA toll free), 321-821-9433 or visit www.AironUSA.com.

ABOUT AIRON CORPORATION

Airon develops, manufactures and distributes pneumatic respiratory products for emergency medicine and critical care. Airon's core products are ventilators and CPAP systems for use in hospitals, MRI, air transport, and emergency medical systems. These products are based upon our patented "High Tech Pneumatics" which has advanced clinical care standards using oxygen-powered devices. The pNeuton ventilators have been on the USA market since 2003 and are now sold in 42 countries around the world. All of Airon's products have FDA marketing clearance and CE mark. Airon is an ISO 9001 and ISO 13485 certified quality manufacturer. For more information on Airon Corporation, please visit www.AironUSA.com.

###

Airon Corporation
751 North Drive, Unit 6, Melbourne, FL 32934 USA
tel (321) 821-9433 fax (321) 821-9443