

pNeuton

Pneumatic Transport Ventilation

Preparing for a disaster?



pNeuton model A for adults / pediatric



pNeuton Mini for neonates / infants

pNeuton is the ideal choice in ventilators.....

- ✓ **Meets USA clinical requirements for disaster preparation** – critically ill patients will require the clinical performance suggested by the CDC including volume or pressure targeted ventilation with PEEP / CPAP and alarms.
- ✓ **Provides invasive or non-invasive ventilation** – a variety of respiratory support options may be required. Not all flu / disaster patients need intubation, face mask CPAP may be appropriate for ventilatory support.
- ✓ **Is easy to use** – designed for inexperienced caregivers. Can even be used by staff that are being remotely supervised by physicians / licensed healthcare providers. Standard tidal volume and rate controls with integrated PEEP / CPAP allow for immediate set-up and extended patient support. Training videos are available.
- ✓ **Works without electricity or batteries** – fully pneumatic system does not require electricity to operate and will function indefinitely if the power goes out. This also allows for unlimited storage without maintenance. Save over \$1,000 per year in storage and maintenance cost compared to battery based ventilators.
- ✓ **Made in the USA** – produced from 95% American made parts in Melbourne, Florida. Service and parts are always available.
- ✓ **Remote alarm compatible** – allows the ventilator to be attached to standard remote alarm nurse call systems. In large isolation rooms with multiple patients, any patient disconnection alarm can be sent to a central nurse station.
- ✓ **MRI compatible** – safe to use within the MRI suite or other special procedures.



Airon

Airon Corporation
751 North Drive, Unit 6, Melbourne, FL 32934 USA
tel +1-321-821-9433 fax +1-321-821-9443
CD-A-041 Rev B