



Factory Rebuild

All medical devices have a defined life cycle for quality, safety and continued effectiveness in meeting performance specifications. With established Preventative Maintenance procedure performed every 2 years, **pNeuton** Ventilators and **MACS** CPAP Systems have an expected maximum 10 year life. For the safety of continued product use beyond its expected life, Airon Corporation recommends a Factory Rebuild manufacturing process. The rebuild is designed to ensure proper performance of the device; your ventilator / CPAP system will look like new and come with a one year warranty.

A complete rebuild includes:

- Replacement of enclosure, pressure gauge and all control knobs so the device will look brand new
- Installation of the PM kit which contains various components including tubing, filters, check valves and reed alarm
- Replacement of all key internal components, ensuring all parts are up to the latest specifications
- Complete inspection of the device
- Calibration and testing to meet factory specifications and ensure your warranty
- Keep the original Serial Number for ease of tracking

All rebuild processes are performed by Airon Corporation at our Melbourne, Florida facility.

Purchase orders for the Rebuild must be issued in advance of the workmanship. All orders must be placed on a direct basis to Airon, not through Airon distributors. Upon receipt of a purchase order, a Returned Goods Authorization Number (RGA) will be issued for the unit to be shipped to Airon. No device will be accepted without a properly issued RGA number. Please contact Technical Services at **1-888-448-1238** or technicalsupport@AironUSA.com to schedule a Rebuild and obtain an RGA.

Prices*:	pNeuton S Rebuild	\$ 3,900.00 – Part number 20002R
	pNeuton A Rebuild	\$ 4,900.00 – Part number 20052R
	MACS CPAP System Rebuild	\$ 1,250.00 – Part number 20201R

Shipping is at the cost of the customer. Upon receipt of the device, we average a 15 business day turn around for the Rebuild. Please verify the estimated service time when requesting a RGA with your purchase order. Loaners are available if desired, at an additional cost.

pNeuton S Loaner	\$ 110.00 – Part number 70070
pNeuton A Loaner	\$ 175.00 – Part number 70071
MACS CPAP System Loaner	\$ 90.00 – Part number 70072

Terms and Conditions are subject to change at any time without notice.

* All pricing is listed for USA only.

Airon Corporation
751 North Drive, Unit 6, Melbourne, FL 32934 USA
tel (321) 821-9433 fax (321) 821-9443